


Quartz crystal microbalance technology as a new platform for real-time monitoring in biology, medicine and engineering

Dr. Paula Braun, Jin Zhang &
Dr. Frank K. Gehring, 3T-analytik, GmbH & Co.KG,
Gartenstr. 100, DE-78532 Tuttlingen, Germany
web: www.3t-analytik.de
email: info@3t-analytik.de

Quartz crystal microbalance (QCM-D) technology is a surface-sensitive technique of (bio)-layers on a surface with regard to adsorption/desorption events, molecular interactions and structural properties. Here, three examples illustrate the versatility of the QCM-D technology in the analysis complex dynamic processes.

- 1) Build-up of polyelectrolyte multi-layers
- 2) Tissue engineering of bio-mimetic bone matrix construct
- 3) Formation and removal of bacterial biofilm layer



Quartz crystal microbalance measuring principle



- At the heart of the QCM-D technology is the biosensor, i.e. a quartz oscillator (Fig. 1).
- Changes in the QCM-D parameters, frequency and damping of the oscillation, correlate with binding to and/or with changes in the structural properties at the surface layer (1, 2).

Three distinct typical cases illustrate the QCM-D measuring principle (Fig. 2):

- Formation of rigid homogeneous layers leads to a change in frequency, with an unchanged damping change in mass but constant viscoelasticity (2A).
- A change in the viscosity of pure liquids, for example, water compared with glycerol, results in a decrease in the frequency while the same amount is increased in the damping (2B).
- The addition of large 'soft' molecules, such as immune cells or bacterial cells leads to a frequency and damping change in opposite directions, as yet, with different amplitudes (2C).
- In particular, the damping parameter permits the detection of events or transformations in 'soft' cellular layers (3) (2C).




Fig. 2 QCM-D experiment on biological samples. Idealized signal shape of frequency (blue line) and damping (yellow line) during three distinct cases. Rigid globular molecules (red circles). Bacterial or immune cells (yellow ovals). See text for further details.

1) Monitoring of polyelectrolyte multilayer build-up

Distinct PM growth and layer thickness




Fig. 3. PM build-up monitored by QCM-D. Progression of the frequency (red line) and the damping (blue line) signals of linearly growing (PEI-(PSS/PAH)₁₀, top) and exponentially growing (PEI-(PGA/PAH)₈, bottom) PMs (A). The arrows point to polyelectrolyte (PE) addition (PEI black, PSS red, PAH blue). PSS/PAH and PGA/PAH film thickness on sensor surface as a function of PE addition (B). Thickness is determined by applying the Sauerbrey (grey line) or the viscoelastic module (black line) (1,2). PEI, polyethylenimine; PSS, polystyrene sulfonate; PAH, polyallylamine; PGA, poly-L-glutamic acid.

2) Monitoring biomimetic bone cell-matrix synthesis

Bio-activation of PLC nanofibers

- Bioactive implants intended for rapid, and durable bone tissue regeneration are generated based on nanofibrous 3D-scaffolds of bio-resorbable poly ε-caprolactone (PCL) mimicking the architecture of bone matrix (5).
- The frequency and damping shifts reflect the gradual stepwise coating of the PCL-nanofiber with linker peptides, poly-lysine (pLys) and – glutamate (pGlu) and at the final coating stage, replenishment of the nanostructures' cisterns with fibroblasts growth factor for bio-activation (Fig. 4, red arrow) (5,6).





Fig. 4 Tissue engineering with QCM-D. (A) High resolution imaging of PCL polymer, coated with the osteogenetic growth factor bone morphogenetic protein 2 (BMP-2). (B) QCM-D frequency signal progression during the stepwise coating of the nanofiber scaffold with, pLys (blue arrows), pGlu (yellow arrows) and the final bio-activation by addition of fibroblast growth factor, FGF2 (red arrow) (6).

3) Real time monitoring of biofilm development

Pseudomonas aeruginosa surface adhesion, biofilm formation and removal



- A portfolio of effective measures is much needed to contain the medically significant biofilms of e.g. multi-drug resistant *Pseudomonas aeruginosa* (*P. aeruginosa*) (7).

- Key information about the function, community diversity, assembly and disassembly, as well as the resistance to disinfectants is obtained by QCM-D monitoring of biofilm maturation (8,9; Fig. 5).

Fig. 5 Biofilm formation monitored with QCM-D. Electron microscopy image of *P. aeruginosa*, http://phil.cdc.gov/phil/details_linked.asp?pid=16876 (A). Adhesion, growth and removal of *P. aeruginosa* monitored by QCM-D-technology (B). Frequency (blue) and damping signal (red) progression. The inset shows the process of adhesion after bacteria inlet (black arrow) during first 60 min. Arrows mark the addition of bacteria (I), growth medium (II) and detergent (SDS) (III). Model scheme of biofilm formation (C): I) Bacteria adhesion, II) growth of the bacterial biofilm and increase in the hydro-viscous state due to emergence of the extra polymeric substances; III) removal of bacterial film from the sensor.

1. Sauerbrey G. 1959 Verwendung von Schwingquarzen zur Wägung dünner Schichten und zur Mikrowägung. *Z Phys*, 155: 206–222. 2. Kanazawa K, Gordon JG. 1985 The oscillation frequency of a quartz resonator in contact with a liquid. *Anal Chim Acta*, 175: 99–105. 3. Gehring, FK. 2006 Schwingquarzensensorik in Flüssigkeiten – Entwicklung eines Blutanalysegerätes, Cuvillier Verlag, Göttingen, ISBN 9783865378729. 4. Grieshaber J, Voros T, Zambelli V, Ball P, Schaaf J-C, Voegel, Boulmedais F. 2008 *Langmuir* 24: 13668–13676. 5. Eap S, Richert L, Lemoine S, Kalaskar D, Demoustier-Champagne S, Atmani H, Mély Y, Fioretti F, Schlatter G, Kuhn L, Ladam G, Benkirane-Jessel N. 2014 Osteogenic properties of electrospun nanofibrous PCL scaffolds equipped with chitosan-based nanoreservoirs of growth factors. *Macromol Biosci*. 14: 45–55. 7. Waslathanthri DP, Kuhn L, Rusling JF. 2013 Multilayer Bio-Films for bone tissue regeneration. University of Connecticut Health Center, Storrs CT. *Appl. Note: www.3t-analytik.de/sites/default/files/Multilayer_Bio-films_for_Bone_Tissue_Regeneration30032014.pdf*. 8. Boucher HW, Talbot GH, Bradley JS et al. 2009 Bad bugs, no drugs: no ESKAPE! An update from the Infectious Diseases Society of America. *Clin Infect Dis* 48:1–12. 8. Olsson AL, Mitzel MR, Tufenkji N. 2015 QCM-D for non-destructive real-time assessment of *Pseudomonas aeruginosa* biofilm attachment to the substratum during biofilm growth. *Colloids Surf B Biointerfaces*, 136: 928–934. 9. Sismaet HJ, Abadian PN, Goluch ED. 2014 Monitoring Bacterial Biofilm Growth and Removal. Department of Chemical Engineering, Northeastern University, Boston, MA, USA. *Appl. Note: www.3t-analytik.de/sites/default/files/AN2015_Bacterial_Biofilm_Growth_and_Removal_.pdf*.